



# Micro<sup>+</sup> Smokerlyzer<sup>®</sup>

Helping people to stop smoking, one breath at a time.



One CO monitor for all  
your testing needs

**FREE**  
COdata<sup>+</sup><sup>™</sup>  
software

Ideal for in-depth stop smoking  
advice delivered by:

- Stop smoking services
- General practice
- Schools
- Pharmacies
- Health visitors
- Midwives
- Researchers



Automatic %FCOHb  
conversion



Create and save  
patient details



For ambient levels  
of CO

**Smokerlyzer<sup>®</sup>**  
[www.bedfont.com](http://www.bedfont.com)

## Features and benefits

- Brand new easy to use interface
- Create and store patient details
- SteriTouch® technology for improve infection control
- Large touch screen
- Automatic conversion to %FCOHb
- Familiar green, amber and red traffic light system
- Ambient levels of CO

## Technical specification

Concentration range	0-500ppm
Display	Full colour touchscreen
Detection principle	Electrochemical sensor
Repeatability	<±5%
Accuracy	2ppm/±5%*
Power	3 x AA (LR6 or equivalent) - up to 1000 minutes 1 x CR2032 Lithium coin cell
T <sub>90</sub> response time	<30 seconds
Operating temperature	15-40°C
Storage/transport temperature	0-50°C
Operating/storage/transport pressure	Atmospheric ±10%
Operating humidity	15-90% non-condensing
Storage/transport humidity	0-95%
Sensor operating life	5 years
Sensor sensitivity	1ppm
Sensor drift	<5% per annum
Dimensions	Approx. 37 x 77 x 140mm
Weight	Approx. 215g (including batteries)
Materials	Case: polycarbonate/ABS blend with SteriTouch® anti-microbial additive D-piece™: polypropylene SteriBreath™: polypropylene
H <sub>2</sub> cross interference	≤6%

\*Whichever is greater



## British technology

**Bedfont® Scientific Ltd**

Station Road, Harrietsham, Maidstone

Kent, ME17 1JA, England

Tel: +44 (0)1622 851122, Fax: +44 (0)1622 854860

Email: [ask@bedfont.com](mailto:ask@bedfont.com)

[www.bedfont.com](http://www.bedfont.com)

© Bedfont® Scientific Limited 2017

Issue 8- April 2017- Part No: MKT375

Bedfont® Scientific Limited reserve the right to change or update this literature without prior notice.

Registered in: England and Wales. Registered No: 1289798



**bedfont®**  
est. 1976